




TMJ YANG'S SYSTEM

Product Catalog of Customized
TMJ-Cranio-Jaw Prosthesis



PRODUCT CATALOG

Customized Total TMJ Prosthesis	01
Customized TMJ-Skull Base Combined Prosthesis	05
Customized Skull Base-Articular Fossa Prosthesis	09
Customized TMJ-Mandibular Combined Prosthesis	13
Standard Total TMJ Prosthesis	17


Customized Total TMJ Prosthesis




Classic Total TMJ Prosthesis





Total TMJ Prosthesis: Mandibular Positioning





Total TMJ Prosthesis: Metal/Non-metal Fossa Prosthesis






Preoperative Examination





Surgical Planning




Surgical Procedure




Surgical Outcome




Customized TMJ-Skull Base Combined Prosthesis





TMJ-Skull Base Combined Prosthesis: Peri-fixed





TMJ-Skull Base Combined Prosthesis: Fenestrated Bone Plate-fixed






Preoperative Examination





Surgical Planning



Surgical Procedure





Surgical Outcome





Customized Skull Base-Fossa Prosthesis





Skull Base-Fossa Prosthesis






Preoperative Examination




Surgical Planning




Surgical Procedure




Surgical Outcome




Customized TMJ-Mandibular Combined Prosthesis





TMJ-Mandibular Combined Prosthesis: Wide-body





TMJ-Mandibular Combined Prosthesis:
Mandibular Body Augmentation





Preoperative Examination





Surgical Planning



Surgical Procedure




Surgical Outcome



Standard Total TMJ Prosthesis


Standard TMJ Prosthesis

Model	l/mm	h/mm
GFGJL26H45	26	45
GFGJL26H50	26	50
GFGJL26H55	26	55




Standard TMJ Prosthesis

Model	l/mm	h/mm
GFGJL30H45	30	45
GFGJL30H50	30	50
GFGJL30H55	30	55



Bone Screw, Cross-Head

Model	d ₁ /mm	l/mm
2.0 MP Self-tapping screw	2.0	3~20 (1 mm increments per specification)
2.1 MP Self-tapping screw	2.1	
2.4 MP Self-tapping screw	2.4	
2.7 MP Self-tapping screw	2.7	
2.9 MP Self-tapping screw	2.9	





Gaofeng Medical WeChat Official Account